

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Rabbit Serum IgG (H+L), DEAE Purified

Catalog Number: 70469 (2002)

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Form: Liquid

Purity: Greater than 98% pure

Buffer: 15 mM potassium phosphate buffer, 0.85% sodium chloride, pH=7.4

Preservative: 0.10% sodium azide
This product is filtered through a 0.2 micron membrane.

Storage: 2-8°C