

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: *Toxoplasma gondii* IgG Antigen

Catalog Number: 70494

Lot Number:

Protein Concentration:

Quantity:

Strain: RH

Produced In: Mice (CD-1)

Product Form: In 20 mM Phosphate buffer, pH=7.40, containing 0.85% NaCl, 1 mM AEBSF, and 10 KIU Aprotinin/ml.

Infectivity: Negative. Gamma Irradiated (30 K rad)

Storage Condition: -70°C or below

Suggested Use: For research use only

Caution: This material should be handled as if it were capable of transmitting infection. Please use good laboratory practice techniques.