

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Cytomegalovirus *IgG Antigen*

Catalog Number: 70537

Lot Number:

Protein Concentration:

Quantity:

Strain: AD169

Produced In: MRC-5 Cells

Product Form: In 20 mM Phosphate Buffer Containing 140 mM NaCl, pH = 7.4

Infectivity: Negative. Gamma Irradiated (30 K rad)

Detergents: < 0.1% CHAPS

Storage Condition: -70°C or below

Suggested Use: For Research Use Only

Caution: This material should be handled as if it were capable of transmitting infection. Please use good laboratory practice techniques.