

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Herpes Simplex Virus, Type 1, IgG Antigen

**Catalog Number:** 70539

**Lot Number:**

**Protein Concentration:**

**Quantity:**

**Strain:** Recombinant Antigen

**Produced In:** E. Coli

**Product Form:** In 20 mM phosphate buffer, 500 mM NaCl, pH = 8.0

**Infectivity:** Negative

**Storage Condition:** -70°C or below

**Suggested Use:** For Research Use Only

**Caution:** This material has been engineered and expressed without the use of whole virus. It contains no HSV-1 virus. However, as with all biological materials, it should be handle as potentially infections.