

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Mouse Monoclonal Anti-CA15-3

Catalog Number: 70563

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Subclass: IgG2b

Product Form: Protein A column chromatography purified. In 0.015 M KPO_4 buffer, pH = 7.40, containing 0.85% NaCl and 0.05% NaN_3

Specificity: CA 15-3 antigen

Storage: At 2-8° C. DO NOT FREEZE. Precipitation may occur upon freezing.

Suggested Use: For labeling use (detection)

Pairing Information: This MoAb can be used as detection antibody when paired with CA15-3 MoAb 70642 or 70650 for solid phase coating (capture).