

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Mouse Monoclonal Anti-Human GI Cancer Antigen
(CA 19-9)-HRP-Conjugate

Catalog Number: 70564-HRP

Lot Number:

Quantity: _____ml

Antibody Concentration (A_{280nm}): _____mg/ml

Product Form: In 0.05 M KPO_4 buffer, 0.85% NaCl, pH=7.20.
Containing protein stabilizers and preservatives

Storage: At 2-8°C. **Do not freeze.**

Dilution Factor: Greater than 1:1,000 in the human CA 19-9 ELISA assay