

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Mouse Monoclonal Anti-Human CA 19-9 Tumor Marker

Catalog Number: 70564

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Subclass: IgG₁

Product Form: Protein A column chromatography purified. In 0.015 M KPO₄ buffer, pH = 7.40, containing 0.85% NaCl and 0.05% NaN₃

Specificity: CA 19-9 antigen

Storage: At 2-8° C. DO NOT FREEZE. Precipitation may occur upon freezing.

Suggested Use: For labeling (detection) use

Pairing Information: This antibody pairs with CA 19-9 MoAb 70576 in the human CA 19-9 sandwich immunoassay.