

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Mouse Monoclonal Anti-Human CA 19-9 Tumor Marker

**Catalog Number:** 70576

**Lot Number:**

**Quantity:** \_\_\_\_\_ mg (= \_\_\_\_\_ ml)

**Antibody Concentration ( $A_{280nm}$ ):** \_\_\_\_\_ mg/ml

**Subclass:** IgG<sub>1</sub>

**Product Form:** Protein A column chromatography purified. In 0.015 M KPO<sub>4</sub> buffer, pH = 7.40, containing 0.85% NaCl and 0.05% NaN<sub>3</sub>

**Specificity:** CA 19-9 antigen

**Storage:** At 2-8° C. DO NOT FREEZE. Precipitation may occur upon freezing.

**Suggested Use:** For solid phase coating (capture)

**Pairing Information:** This antibody pairs with CA 19-9 MoAb (P/N: 70564) in the human CA 19-9 sandwich immunoassay.