

BIOCHECK, INC.
323 Vintage Park Drive
Foster City, CA 94404
Tel: 650-573-1968 Fax: 650-573-1969

PRODUCT SPECIFICATIONS

Product Name: Rabbit Polyclonal Anti-Testosterone, Neat Serum

Catalog Number: 70616

Antigen Source: Testosterone-19-HSA

Lot Number:

Quantity: _____ ml (Neat Serum)

Product Form: Rabbit Anti-Testosterone, Neat Serum, 0.1% NaN₃

Storage: -20°C

Cross-Reactivity:

Testosterone	100%
Dihydro Testosterone	3.5%
Danazol	0.004%
Nortestosterone	<17.0%
Sex Hormone Binding Globulin	0%
Other Steroids	<0.01%

Suggested Use: For Testosterone Competitive Immunoassay.