

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Mouse Monoclonal Anti-Testosterone

Immunogen: Testosterone-3-CMO-BSA

Catalog Number: 70619

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: DEAE column chromatography purified immunoglobulin fractions in 0.05 M PBS, pH = 7.40, containing 0.1% NaN₃.

Affinity Constant: 2×10^9 l/mol.

Subclass: IgG_{2a}

Storage: At 2-8° C. DO NOT FREEZE. Precipitation may occur upon freezing.

Specificity: Testosterone 100%
Very low cross reactivity to other steroids.