

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Goat Polyclonal Anti-Human-Whole-TSH, Affinity Purified

Catalog Number: 70646

Lot Number:

Quantity: _____ mg (= _____ ml)

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: TSH antigen affinity columns purified immunoglobulins in 0.015 M KPO_4 buffer, pH=7.40, containing 0.85% NaCl and 0.1% NaN_3

Storage: At 2-8° C. DO NOT FREEZE. Precipitation may occur upon freezing.

Specificity: Human Whole-TSH, alpha-TSH, and beta-TSH

Suggested Use: For enzyme conjugation (detection) or solid phase coating (capture).

Pairing Information: This antibody pairs with TSH MoAb (P/N: 70204) in the human TSH sandwich immunoassay.