

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Rabbit Polyclonal Anti-Human CRP, Affinity Purified.

**Code:** 70648

**Lot Number:**

**Quantity:** \_\_\_\_\_ mg (= \_\_\_\_\_ ml )

**Antibody Concentration ( $A_{280nm}$ ):** \_\_\_\_\_ mg/ml

**Product Form:** CRP antigen affinity column chromatography purified antibodies in 15 mM potassium phosphate buffer, pH = 7.40, containing 0.85% NaCl and 0.1% NaN<sub>3</sub>.

**Storage:** 2-8° C

**Specificity:** This antibody reacts with human CRP.

**Suggested Use:** For enzyme conjugation (detection).

**Pairing Information:** This antibody pairs with goat anti-CRP (P/N: 70243) in the human CRP sandwich immunoassay.