

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name:	Mouse Monoclonal Anti-CA15-3
Catalog Number:	70650
Lot Number:	
Quantity:	_____ mg (= _____ ml)
Antibody Concentration (A_{280nm}):	_____ mg/ml
Subclass:	IgG1
Product Form:	Protein A column chromatography purified. In 0.015 M KPO_4 buffer, pH = 7.40, containing 0.85% NaCl and 0.05% NaN_3
Specificity:	This monoclonal antibody is specific to carbohydrate on MUC-1.
Storage:	At 2-8° C. <u>DO NOT FREEZE</u> . Precipitation may occur upon freezing.
Suggested Use:	For solid phase coating at pH=6.5 (capture)
Pairing Information:	This antibody can be used as solid phase coating (capture) when paired with CA15-3 MoAb 70563 as detection antibody.
Important Note:	For solid phase coating, please use 0.2M $NaPO_4$ buffer, pH=6.50. This antibody is not stable at pH= 9.60.