

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name:	Mouse Monoclonal Anti-Human Hemoglobin
Catalog Number:	70658
Lot Number:	
Quantity:	_____ mg (= _____ ml)
Antibody Concentration (A_{280nm}):	_____ mg/ml
Subclass:	IgG ₁
Product Form:	Protein A column chromatography purified immunoglobulin fractions in 0.015M KPO ₄ buffer, pH = 7.40, containing 0.85% NaCl and 0.1% NaN ₃ .
Storage:	At 2-8° C. <u>DO NOT FREEZE</u> . Precipitation may occur upon freezing.
Specificity:	Human Hemoglobin
Suggested Use:	For enzyme, colloidal gold or latex conjugation in the human hemoglobin sandwich immunoassay.
Pairing Information:	This antibody pairs with the following antibody in the human hemoglobin sandwich Immunoassay. <ul style="list-style-type: none">• Mouse monoclonal anti-human hemoglobin P/N:70657 (capture).