

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name:	Mouse Monoclonal Anti-Human Lipoprotein-Associated Phospholipase A2 (Lp-PLA2)
Catalog Number:	70745
Lot Number:	
Quantity:	
Antibody Concentration (A_{280nm}):	
Subclass:	IgG ₁
Product Form:	PBS, pH 7.4 with 0.09% sodium azide
Storage:	At 2-8°C. <u>Do Not Freeze.</u> Precipitation may occur upon freezing.
Suggested Use:	For solid phase coating (capture) in the human Lp-PLA2 immunoassay.
Pairing Information:	This antibody pairs with Mouse Monoclonal Anti-Human Lp-PLA2 Part No. 70746 or Mouse Monoclonal Anti-Human Lp-PLA2 Part No. 70747