

BIOCHECK, INC.

425 Eccles Avenue
South San Francisco, CA 94080
Tel: (650) 573-1968 E-mail: info@biocheckinc.com

PRODUCT SPECIFICATIONS

Product Name:	Mouse Monoclonal Anti-Human Lipoprotein-Associated Phospholipase A2 (Lp-PLA2)
Catalog Number:	70745
Lot Number:	
Quantity:	
Antibody Concentration (A_{280nm}):	
Subclass:	IgG ₁
Product Form:	PBS, pH 7.4 with 0.09% sodium azide
Storage:	At 2-8°C. <u>Do Not Freeze.</u> Precipitation may occur upon freezing.
Suggested Use:	For solid phase coating (capture) in the human PD-1 immunoassay.
Pairing Information:	This antibody pairs with Mouse Monoclonal Anti-Human Lp-PLA2 Part No. 70746 or Mouse Monoclonal Anti-Human Lp-PLA2 Part No. 70747