

# BIOCHECK, INC.

425 Eccles Avenue  
South San Francisco, CA 94080  
Tel: (650) 573-1968 E-mail: info@biocheckinc.com

---

## PRODUCT SPECIFICATIONS

<b>Product Name:</b>	<b>Mouse Monoclonal Anti-Human PD-1</b>
<b>Catalog Number:</b>	70749
<b>Quantity:</b>	0.1 mg, 0.5 mg or 1.0 mg
<b>Subclass:</b>	IgG <sub>1</sub>
<b>Product Form:</b>	Lyophilized (Carrier-Free)
<b>Storage:</b>	Store at 2-8°C.
<b>Suggested Use:</b>	For solid phase coating ( <b>capture</b> ) in the human PD-1 immunoassay.
<b>Pairing Information:</b>	This antibody pairs with Mouse Monoclonal Anti-PD-1 Part No. 70750.