

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name:	Mouse Monoclonal Anti-Human PD-1
Catalog Number:	70750
Quantity:	0.1 mg, 0.5 mg, or 1.0 mg
Subclass:	IgG2a
Product Form:	Lyophilized (Carrier-Free)
Storage:	Store at 2-8°C.
Suggested Use:	For enzyme conjugation (detection) in the PD-1 immunoassay.
Pairing Information:	This antibody pairs with Mouse Monoclonal Anti-PD-1 Part No. 70748 or Mouse Monoclonal Anti-PD-1 Part No. 70749.