

BIOCHECK, INC.

425 Eccles Avenue
South San Francisco, CA 94080
Tel: (650) 573-1968 E-mail: info@biocheckinc.com

PRODUCT SPECIFICATIONS

Product Name:	Mouse Monoclonal Anti-Human PD-L1
Catalog Number:	70751
Quantity:	0.1 mg, 0.5 mg, or 1.0 mg
Subclass:	IgG ₁
Product Form:	Lyophilized (Carrier-Free)
Storage:	Store at 2-8°C.
Suggested Use:	For solid phase coating (capture) in the human PD-L1 immunoassay.
Pairing Information:	This antibody pairs with Mouse Monoclonal Anti-PD-L1 Part No. 70753.