

BIOCHECK, INC.

425 Eccles Avenue
South San Francisco, CA 94080
Tel: (650) 573-1968 E-mail: info@biocheckinc.com

PRODUCT SPECIFICATIONS

| | |
|-----------------------------|--|
| Product Name: | Mouse Monoclonal Anti-Human PD-L1 |
| Catalog Number: | 70752 |
| Quantity: | 0.1 mg, 0.5 mg, or 1.0 mg |
| Subclass: | IgG ₁ |
| Product Form: | Lyophilized (Carrier-Free) |
| Storage: | Store at 2-8°C. |
| Suggested Use: | For solid phase coating (capture) in the human PD-L1 immunoassay. |
| Pairing Information: | This antibody pairs with Mouse Monoclonal Anti-PD-L1 Part No. 70753. |