

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Goat Polyclonal Anti-Human Beta-TSH-HRP Conjugate

Catalog Number: 80001 (G-7003-HRP)

Lot Number:

Quantity: _____ ml

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: In 0.05 M KPO_4 buffer, 0.85% NaCl, pH=7.20, containing proteins and preservatives.

Storage: At 2-8° C. DO NOT FREEZE. The conjugate may become inactivated upon freezing.

Pairing Inofrmation: This conjugate pairs with mouse monoclonal anti-Beta-TSH (P/N: 70204) in the **human TSH sandwich immunoassay**.