

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: T3-HRP Conjugate

Catalog Number: 80003

Lot Number:

Quantity: _____ml

Product Form: 1:1 (v/v) (20 mM KPO₄ buffer, 0.85% (w/v) NaCl, pH = 7.2) / glycerol

Storage: -20°C

Suggested Use: For competitive human T3 immunoassay.