

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Mouse Monoclonal Anti-Human Alpha Fetoprotein (AFP) HRP Conjugate

Catalog Number: 80005 (5005-HRP)

Lot Number:

Quantity: _____ ml

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: In 0.05 M KPO_4 buffer, 0.85% NaCl, pH=7.20, containing protein stabilizers and preservatives.

Storage: At 2-8°C. **Do not freeze.**

Dilution Factor: Greater than 1:1,000 in the AFP ELISA assay.

Pairing Information: This conjugate pairs with rabbit polyclonal anti-AFP (P/N: 70269) in the **AFP sandwich immunoassay**.