

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Goat Polyclonal Anti-CEA-HRP Conjugate

Catalog Number: 80006 (3005-HRP)

Lot Number:

Quantity: _____ ml

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: In 0.05M KPO_4 buffer, 0.85% NaCl, pH=7.20, containing proteins and preservatives.

Storage: At 2-8°C. DO NOT FREEZE. The conjugate may become inactivated upon freezing.

Dilution factor: Greater than 1: 1,000 in the CEA ELISA assay.

Pairing Information: This antibody pairs with CEA MoAb (P/N: 70357 and 70370) in **the CEA sandwich immunoassay**.