

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Mouse Monoclonal Anti-Human Breast Cancer Antigen (CA15-3)-HRP-Conjugate

Catalog Number: 80008 (CA15-3 MoAb-HRP)

Lot Number:

Quantity: _____ ml

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: In 0.05 M KPO_4 buffer, 0.85% NaCl, pH=7.20, containing protein stabilizers and preservatives

Storage: At 2-8°C. Do not freeze.

Dilution Factor: Greater than 1:1,000 in the human CA15-3 ELISA assay

Pairing Information: This conjugate pairs with CA15-3 MoAb, P/N: 70642, in the human CA15-3 sandwich immunoassay.