

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Mouse Monoclonal Anti-PSA -HRP Conjugate

**Catalog Number:** 80010 (PSA-103-HRP)

**Lot Number:**

**Amount:** \_\_\_\_\_ ml

**Antibody Concentration ( $A_{280nm}$ ):** \_\_\_\_\_ mg/ml

**Product Form:** In 0.05 M  $KPO_4$  buffer, 0.085% NaCl, pH = 7.20, containing protein stabilizers and preservatives.

**Storage:** At 2-8°C. DO NOT FREEZE. The Conjugate may become inactivated upon freezing.

**Dilution:** Greater than 1:1,000 in the PSA ELISA assay

**Pairing Information:** This conjugate pairs with **goat anti-PSA (P/N: 70415) in the Total PSA sandwich immunoassay**