

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Mouse Monoclonal Anti-Human PSA -HRP Conjugate

Catalog Number: 80010 (PSA-103-HRP)

Lot Number:

Amount: _____ml

Antibody Concentration (A_{280nm}): _____mg/ml

Product Form: In 0.05 M KPO_4 buffer, 0.085% NaCl, pH = 7.20, containing protein stabilizers and preservatives.

Storage: At 2-8°C. DO NOT FREEZE. The Conjugate may become inactivated upon freezing.

Dilution: Greater than 1:1,000 in the human PSA ELISA assay

Pairing Information: This conjugate pairs with goat polyclonal anti-human PSA (P/N: 70415) in the human **Total PSA** sandwich Immunoassay.