

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name:	Goat Polyclonal Anti-Human Prostate Specific Antigen (PSA)- HRP Conjugate
Catalog Number:	80011
Lot Number:	
Quantity:	_____ ml
Antibody Concentration (A_{280nm}):	_____ mg/ml
Product Form:	In 0.05 M KPO ₄ buffer, 0.85% NaCl, pH=7.20, containing proteins and preservatives.
Storage:	At 2-8° C. <u>DO NOT FREEZE</u> . The conjugate may become inactivated upon freezing.
Dilution Factor:	Greater than 1:1,000 in the Free PSA ELISA assay.
Suggested Use:	<ol style="list-style-type: none">1. This conjugate pairs with mouse monoclonal anti-PSA (P/N: 70338) in the <u>Total PSA</u> sandwich immunoassay.2. This conjugates pairs with mouse monoclonal anti-<u>Free PSA</u> (P/N: 70340) in the <u>Free PSA</u> sandwich immunoassay.