

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

|   |  |
|---|--|
| <b>Product Name:</b>                                    | <b>Goat Polyclonal Anti-Human Prostate Specific Antigen (PSA)- HRP Conjugate</b>   |
| <b>Catalog Number:</b>                                  | 80011  |
| <b>Lot Number:</b>                                      |  |
| <b>Quantity:</b>  | _____ ml   |
| <b>Antibody Concentration (<math>A_{280nm}</math>):</b> | _____ mg/ml  |
| <b>Product Form:</b>                                    | In 0.05 M $KPO_4$ buffer, 0.85% NaCl, pH=7.20, containing proteins and preservatives.  |
| <b>Storage:</b>   | At 2-8° C. <u>DO NOT FREEZE</u> . The conjugate may become inactivated upon freezing.  |
| <b>Dilution Factor:</b>                                 | Greater than 1:1,000 in the human Free PSA ELISA assay.  |
| <b>Suggested Use:</b>                                   | <ol style="list-style-type: none"><li>1. This conjugate pairs with mouse monoclonal anti-human PSA (P/N: 70338) in the human <b><u>Total PSA</u></b> sandwich immunoassay.</li><li>2. This conjugates pairs with mouse monoclonal anti-human <b><u>Free PSA</u></b> (P/N: 70340) in the human <b><u>Free PSA</u></b> sandwich immunoassay.</li></ol> |