

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Mouse Monoclonal Anti- β -FSH HRP Conjugate

Catalog Number: 80015 (G-3001-HRP)

Lot Number:

Quantity: _____ ml

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: In 0.05 M KPO_4 buffer, 0.085% NaCl, pH = 7.20, containing proteins and preservatives.

Storage: At 2-8°C. DO NOT FREEZE. The conjugate may become inactivated upon freezing.

Dilution Factor: Greater than 1:1,000 in the human FSH ELISA assay.

Pairing Information: This conjugate pairs with alpha-FSH MoAb (P/N: 70327) In the human FSH sandwich immunoassay.