

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Mouse Monoclonal Anti-Human Prolactin- HRP Conjugate

**Catalog Number:** 80019 (G-5004-HRP)

**Lot Number:**

**Quantity:** \_\_\_\_\_ ml

**Antibody Concentration ( $A_{280nm}$ )** \_\_\_\_\_ mg/ml

**Product Form:** In 0.05 M  $KPO_4$  buffer, 0.85% NaCl, pH = 7.20, containing proteins and preservatives.

**Storage:** At 2-8°C. DO NOT FREEZE. The conjugate may become inactivated upon freezing.

**Dilution Factor:** Greater than 1:1,000 in the Prolactin ELISA assay.

**Pairing Information:** This conjugate pairs with Prolactin MoAb (P/N: 70317) in **the human Prolactin sandwich immunoassay.**