

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Mouse Monoclonal Anti- $\beta$ -hCG HRP Conjugate

**Catalog Number:** 80023 (5008-HRP)

**Lot Number:**

**Quantity:** \_\_\_\_\_ ml

**Antibody Concentration ( $A_{280\text{nm}}$ ):** \_\_\_\_\_ mg/ml

**Product Form:** In 0.05 M  $\text{KPO}_4$  buffer, 0.085% NaCl, pH = 7.20, containing proteins and preservatives.

**Storage:** At 2-8°C. DO NOT FREEZE. The conjugate may become inactivated upon freezing.

**Dilution Factor:** Greater than 1:1,000 in the hCG ELISA assay.

**Pairing Information:** This conjugate pairs with Alpha hCG MoAb (P/N: 70327) in the hCG sandwich immunoassay.