

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Sheep Polyclonal Anti-Human β -2 Microglobulin-HRP Conjugate

Catalog Number: 80027

Lot Number:

Quantity: _____ ml

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: In 0.05 M KPO_4 buffer, 0.85% NaCl, pH=7.20, containing proteins and preservatives.

Storage: At 2-8° C. DO NOT FREEZE. The conjugate may become inactivated upon freezing.

Dilution Factor: Greater than 1:1,000 in the β -2 microglobulin ELISA assay.

Pairing Information: This conjugate pairs with mouse monoclonal anti- β -2 microglobulin (P/N: 70191) in the **human β -2 microglobulin** sandwich immunoassay.