

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Mouse Monoclonal Anti-Human Cardiac Troponin I MoAb (TPC-110)-HRP Conjugate

Catalog Number: 80055 (TPC-110-HRP)

Lot Number:

Quantity: _____ ml

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: In 0.05M KPO_4 buffer, 0.85% NaCl, PH= 7.20, containing proteins and preservatives.

Storage: At 2-8°C. Do Not Freeze. The conjugate may become inactivated upon freezing.

Specificity: Human Cardiac Specific Troponin I.

Dilution Factor: Greater than 1:1,000 in the Troponin I ELISA assay.

Pairing Information: This conjugate pairs with Troponin I MoAb mixtures (TPC-6+TPC-102+TPC-302 2:2:1 w/w/w) in the Troponin I sandwich immunoassay.