

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Goat Polyclonal Anti-Human Myoglobin  
Antibody- HRP Conjugate

**Catalog Number:** 80058

**Lot Number:**

**Quantity:** \_\_\_\_\_ ml

**Antibody Concentration ( $A_{280nm}$ ):** \_\_\_\_\_ mg/ml

**Product Form:** In 0.05 M  $KPO_4$  buffer, 0.85% NaCl, pH=7.20,  
containing proteins and preservatives.

**Storage:** At 2-8° C. DO NOT FREEZE. The conjugate may  
become inactivated upon freezing.

**Dilution Factor:** Greater than 1:4,000 in the Myoglobin ELISA assay.

**Pairing Information:** This conjugate pairs with mouse monoclonal anti-  
myoglobin (P/N: 70131) in the **myoglobin**  
**sandwich immunoassay**.