

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Goat Polyclonal Anti-Human CRP
Antibody- HRP Conjugate

Catalog Number: 80059

Lot Number:

Quantity: _____ ml

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: In 0.05 M KPO_4 buffer, 0.85% NaCl, pH=7.20,
containing proteins and preservatives.

Storage: At 2-8° C. DO NOT FREEZE. The conjugate may
become inactivated upon freezing.

Dilution Factor: Greater than 1:4,000 in the CRP ELISA assay.

Pairing Information: This conjugate pairs with goat polyclonal anti-CRP
(P/N: 70243) or rabbit polyclonal anti-CRP (P/N: 70324)
in the **CRP sandwich immunoassay**.