

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Goat Anti-Fatty Acid Binding Protein (FABP)  
-HRP Conjugate

**Catalog Number:** 80099

**Lot Number:**

**Quantity:** \_\_\_\_\_ ml

**Antibody Concentration ( $A_{280\text{nm}}$ ):** \_\_\_\_\_ mg/ml

**Product Form:** In 0.015 M  $\text{KPO}_4$  Buffer, 0.85 NaCl, pH = 7.20,  
containing proteins and preservatives.

**Storage:** At 2-8° C. DO NOT FREEZE. The conjugate may  
become inactivated upon freezing.

**Pairing Information:** This conjugate pairs with goat anti-fatty acid binding  
protein (P/N: 70375, for plate coating) in the  
**Fatty Acid Binding Protein** sandwich  
immunoassay.