

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Mouse Monoclonal Anti-Myeloperoxidase (MPO)
-HRP Conjugate

Catalog Number: 80108

Lot Number:

Quantity: _____ ml

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: In 0.05 M KPO_4 buffer, 0.85% NaCl, pH=7.20,
containing proteins and preservatives.

Storage: At 2-8° C. DO NOT FREEZE. The conjugate may
become inactivated upon freezing.

Pairing Information: This conjugate pairs with mouse monoclonal anti-MPO
(P/N: 70355, for plate coating) in the
Myeloperoxidase sandwich immunoassay.