

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICAITONS

Product Name: Mouse Monoclonal Anti-Human PSA-HRP Conjugate

Catalog Number: 81017 (PSA-017-HRP)

Lot Number:

Amount: _____ ml

Antibody Concentration (A_{280nm}): _____ mg/ml

Product Form: In 0.05 M KPO_4 buffer, 0.085% NaCl, pH = 7.20, containing protein stabilizers and preservatives.

Storage: At 2-8°C. DO NOT FREEZE. The Conjugate may become inactivated upon freezing.

Dilution: Greater than 1:1,000 in the human PSA/ACT Complex ELISA assay

Pairing Information: This conjugate pairs with mouse monoclonal anti-human PSA/ACT Complex (P/N: 70349) in the human **PSA/ACT Complex** sandwich immunoassay.