

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Human Thyroid Stimulating Hormone (TSH)
Standard Stock Solution

Catalog Number: 90001

Lot Number:

Specific Immunological Activity: _____ mIU/ml

Quantity: _____ ml

Reference Preparation: WHO, 2nd IRP 80/558

Source: Human Pituitary gland

Storage: - 20°C or below.

Product Form: Prepared in Fetal bovine serum containing 0.5%(v/v)
Proclin-300 as preservative.

Reference Assay: BioCheck TSH ELISA Test

Viral Tests: The source materials used for manufacture of this product tested negative for HBsAg, HIV 1&2 and HCV by FDA-approved methods. No methods can completely assure the absence of these agents and other infectious agents may be present.