

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: CA 15-3 Standard Stock Solution

Catalog Number: 90008

Source: Human BTA cell line supernatant

Lot Number:

CA 15-3 Concentration: _____ U/ml

Quantity: _____ ml

Product Form: In PBS buffer with 0.1% NaN₃, pH = 7.4

Storage: -20°C or below.

Viral Tests: The source materials used for manufacture of this product tested negative for HBsAg, HIV 1&2 and HCV by FDA-approved methods. No methods can completely assure the absence of these agents and other infectious agents may be present.