

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Human Free PSA Standard Stock Solution

**Catalog Number:** 90011

**Source:** Human seminal Fluids

**Lot Number:**

**Free PSA Concentration:** \_\_\_\_\_ ng/ml

**Quantity:** \_\_\_\_\_ ml

**Product Form:** In buffer solution containing 1.0% (v/v) ProClin-300 as preservative.

**Storage:** -20°C or below.

**Viral Tests:** The source materials used for manufacture of this product tested negative for HBsAg, HIV 1&2 and HCV by FDA-approved methods. No methods can completely assure the absence of these agents and other infectious agents may be present.