

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Human Chorionic Gonadotropin, Free Beta-Subunit,
Standard Stock Solution

Catalog Number: 90012

Source: Human Pregnancy Urine

Lot Number:

β -hCG Concentration: _____ ng/ml

Quantity: _____ ml

Product Form: In buffer solution containing 1.0% (v/v) Proclin-300
as preservative

Specific Immunological Activity: 1 ng/ml = 0.1 mIU/ml of hCG

Storage: -20°C or below

Viral Tests: The source materials used for manufacture of this product
tested negative for HBsAg, HIV 1&2 and HCV by FDA-
approved methods. No methods can completely assure the
absence of these agents and other infectious agents may be
present.