

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** Human Chorionic Gonadotropin, Free Beta-Subunit, Standard Stock Solution

**Catalog Number:** 90012

**Lot Number:**

**Source:** Human Pregnancy Urine

**$\beta$ -hCG Concentration:** \_\_\_\_\_ ng/ml

**Quantity:** \_\_\_\_\_ ml

**Product Form:** In newborn bovine serum containing 1.0% (v/v) Proclin-300 as preservative

**Specific Immunological Activity:** 1 ng/ml = 0.1 mIU/ml of hCG

**Storage:** -20°C or below