

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: FSH Standard Stock Solution

Catalog Number: 90015

Source: Human Pituitary Glands

Lot Number:

FSH Concentration: _____ mIU/ml (WHO 2nd IRP 78/549)

Quantity: _____ ml

Product Form: In bovine serum containing 1.0% (v/v) ProClin-300 as preservative.

Storage: -20°C or below.

Viral Tests: The source materials used for manufacture of this product tested negative for HBsAg, HIV 1&2 and HCV by FDA-approved methods. No methods can completely assure the absence of these agents and other infectious agents may be present.