

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Human Luteinizing Hormone (LH),
Standard Stock Solution

Catalog Number: 90016

Lot Number:

Source: Human Pituitary gland

Specific Immunological Activity: _____ mIU/ml

Quantity: _____ ml

Reference Preparation: WHO 1st IRP 68/40

Product Form: In bovine serum containing 1.0% (v/v) ProClin-300 as preservative.

Storage: -20°C or below