

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Troponin I Standard Stock Solution
(Troponin I-T-C Complex)

Catalog Number: 90055 (4003)

Lot Number:

Troponin I Concentration: _____ ng/ml

Quantity: _____ ml

Product Form: In human serum containing 1.0% (v/v) Proclin-300 as preservative.

Storage: -20°C or below.

Reference: Assay: BioCheck Troponin I ELISA Test.

Viral Tests: The source materials used for manufacture of this product tested negative for HBsAg, HIV 1&2 and HCV by FDA-approved methods. No methods can completely assure the absence of these agents and other infectious agents may be present.