

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Myoglobin Standard Stock Solution

Catalog Number: 90058

Source: Human Heart

Lot Number:

Myoglobin Concentration: _____ ng/ml

Quantity: _____ ml

Product Form: 6% BSA in 0.015 M PO₄ buffer, 0.85% NaCl, pH = 7.40
containing 1.0% (v/v) ProClin-300 as preservative.

Storage: -20°C or below.

Viral Tests: The source materials used for manufacture of this product tested negative for HBsAg, HIV 1&2 and HCV by FDA-approved methods. No methods can completely assure the absence of these agents and other infectious agents may be present.