

BIOCHECK, INC.

323 Vintage Park Drive
Foster City, CA 94404
Tel: (650) 573-1968 Fax: (650) 573-1969

PRODUCT SPECIFICATIONS

Product Name: Fatty Acid Binding Protein Standard Stock Solution

Catalog Number: 90099

Lot Number:

Quantity: _____ ml

Specific Immunological Activity: _____ mg/ml

Product Form: Human cardiac tissue

Product Form: In 50% Glycerol, 75 mM NaCl, 5 mM sodium phosphate,
0.05% NaN₃, pH = 7.2

Storage: -20°C or below.