

# BIOCHECK, INC.

323 Vintage Park Drive  
Foster City, CA 94404  
Tel: (650) 573-1968 Fax: (650) 573-1969

---

## PRODUCT SPECIFICATIONS

**Product Name:** *H. pylori* Antigen, Unpurified

**Catalog Number:** HP-1001

**Lot Number:**

**Quantity:** \_\_\_\_\_ mg

**Protein Concentration:** 1.50 mg/ml by Pierce BCA protein assay

**Product Form:** Detergent extract of whole *H. pylori* cell sonicates. In 0.015M KPO<sub>4</sub> buffer, 0.85% NaCl, 0.1% NaN<sub>3</sub>, pH = 7.40.

**Specificity:** *H. pylori* antibody reactive antigens

**Storage Condition:** -20°C or below

**Cross-Reactivity:** *H. pylori* specific. The antigens were extracted from a pure *H. pylori* culture derived from one single *H. pylori* colony.

**Suggested Use:** Microtiter Wells EIA or other immunoassays.

**Caution:** The *H. pylori* antigen is an extract from whole *H. pylori* cells. This material should be handled as if it were capable of transmitting infection. Please use good laboratory practice techniques.